


USB V5 license installation and activation

1/ Installation of Phase Manager


Run the setup and choose the installation language:



Skip the welcome screen by clicking on **Next** :



Accept the license agreement then click on **Next** :



Choose a destination folder **on the USB stick** (for example : E:\PhaseManager31) then click on **Next** :



Choose '**Portable installation**' :

Select Components
Which components should be installed?



Select the components you want to install; clear the components you do not want to install. Click Next when you are ready to continue.

Portable installation ▼

Validate the Start Menu Folder screen by clicking on **Next** :

Select Start Menu Folder
Where should Setup place the program's shortcuts?



 Setup will create the program's shortcuts in the following Start Menu folder.

To continue, click Next. If you would like to select a different folder, click Browse.

Projet Lineaire

< Back **Next >** Cancel

Choose and validate the icons creation screen by clicking **Next** :

Select the additional tasks you would like Setup to perform while installing Phase Manager4 4.0, then click Next.

Additional icons:


☒ Create a desktop icon

☐ Create a Quick Launch icon

On the following screen click on **Install** then once installation is complete, click on **Finish**

2/ USB license configuration


Once Phase Manager 5 is installed, double-click on its icon to launch it.
At first launch, the license configurator will display:



To activate your license, you will need to fill in the following fields:

Module

Choose **USB** :



Product ID

Fill in the Product ID that's been sent to you by Projet Linéaire:

- licence pattern : LPMU-00**NNN**0500abcDEF (**NNN** is your license number, and **5** the version)

- only paste the right part of your Product-ID, after the '-' :

Module	USB	
Product ID	LPMU	00NNN0400abcDEF

Activation key

Fill in the Serial number sent to you by Projet Linéaire.


- pattern : VccQX-37eJG-AZlce-EYk7O-temCP-8DRFu

Module	USB	
Product ID	LPMU	00NNN0400abcDEF
Activation...	VccQX-37eJG-AZlce-EYk7O-temCP-8DRFu	

Licence

Choose “Create new licence” :

Create new licence ☒	Choose a licence ☐
---	--

By clicking on the  button choose a folder **on your USB stick** (for ex. : E:\licence-phasemanager):


Module	USB	
Product ID	LPMU	00NNN0400abcDEF
Activation...	VccQX-37eJG-AZlce-EYk7O-temCP-8DRFu	
Create new licence ☒ Choose a licence ☐		
Browse	E:\License-PM	

3/ License activation

Click now on the  button to activate your licence.

Beware : an Internet connection is mandatory to activate the license !

To activate the license, switch to English then choose the option below and click on **Next** :



The screenshot shows the 'Activation of PhaseManager' window. At the top, it says 'Welcome to PhaseManager'. There are two radio button options. The first option, 'I have a Serial Number and I want to activate PhaseManager', is selected. Below it, a text box prompts the user to 'Enter the Serial Number to activate PhaseManager.'. The second option is 'I want to evaluate PhaseManager (there are 14 day(s) left to evaluate)'. Below this option, it states 'You can evaluate PhaseManager for another 14 day(s). After that you have to activate PhaseManager to further use the product.'. At the bottom left, there is a language dropdown menu set to 'English'. At the bottom right, there are two buttons: 'Next' and 'End'.

- Fill-in all required data then click on Next:
- **The Serial number is the Activation Key that's been used at step 2 in the license configurator**

Activation of PhaseManager

Please enter the activation data

Please enter all required information for the activation process.

Enter Serial Number

JBc F8

[Where do you find your Serial Number?](#)

Enter Name / Company / Organisation

SSA/PL

Enter your E-mail Address


ssauteur@projetlineaire.com

! required


English

Back Next End

- Check that you ARE connected to INTERNET then click on Next:



- If the following message is displayed, that means that Phase Manager did activate right. Click on Next to run Phase Manager:



- If you need to, feel free to contact our customer support at support@projetlineaire.com